

A photograph of a hospital bed with a white sheet and a blue CarePad sensor mat placed under the mattress at the side of the bed. The bed frame is light-colored wood, and the metal rails are visible.

EARLY ALERT IN CASE OF BED ESCAPE

CarePad

THE INVISIBLE ALARM SYSTEM IN THE NURSING BED

APPLICATION

CarePad is a bed escape alarm system and dementia wanderers protection that allows monitoring of „night walkers“ and „wanderers“. The CarePad is placed under the mattress at the side of the bed in the exit zone. When the patient wants to leave the bed, placing the legs on the bed the sensor triggers. Caregivers and/or attendants are alerted.

TECHNICAL DATA

Dimensions: 920 x 150 x 25 mm
Weight: 800 g
Installation: Plug & Play
Connection: compatible with existing alarm and emergency call infrastructure

Areas of application

Retirement and nursing homes, hospitals, clinics, psychiatric wards as well as assisted-living group facilities.

Installation

Unobtrusive placement under the mattress.
All common call systems with protective extra-low voltage or with external power supply

Features

Reliable alarm transmitter with simple operation
Mobile application
No power supply necessary
Robust and easy-to-clean material
Practical housing for radio transmitter

Advantages for nursing staff

Avoidance of falls due to rapid assistance
Reduces night-time checks
No false alarms due to placing the hand or arm, Minimum pressure required.
Unlimited mobility of the caregiver in front of the bed

Advantages for residents or patients

Safety due to invisible protective device
No tripping hazard due to placement under the mattress
Improvement of the quality of life

YOUR RELIABLE PARTNER

Health is a person's most important resource. For this reason, we have set ourselves the task of maintaining it as best as possible and to making an important contribution to recovery. meta-care is serving more than 200 partners and resellers worldwide.



Your CarePad sales representative

Your CarePad distributor

meta-care GmbH

Akazienhofstraße 20
9020 Klagenfurt Austria
Email: office@meta-care.eu

Tel: +43 5 1721 0
Web: www.meta-care.eu